


ENGLISH



HIOSSEN

OSSTEM[®]

IMPLANT

Qualität schafft Vertrauen

Crestal Approach Sinus-KIT

CAS-KIT

CAS-KIT

(Crestal Approach - Sinus KIT)

Contents

- 1 Introduction
- 2 CAS-KIT Feature
- 3 CAS-Drill Specifications & Performance
- 4 Components
- 5 Clinical Indications & Case Study
- 6 Surgical Procedures
- 7 Others




● Introduction ↙

- Hiossen's Crestal Approach Sinus KIT (CAS-KIT) is specifically designed to easily and safely lift the membrane in the maxillary sinus from a crestal approach.

The key component of the CAS-KIT is the CAS-Drill. The unique design of the CAS-Drill enhances convenience and safety of maxillary sinus surgery by; safely lifting the membrane while drilling, precision cutting, flexible cutting speed from low to high speed (800rpm), formation of conical shaped bone chip, generation of bone particles, smooth & stable insertion, easy path correction and septum surgery.

● FEATURES of CAS-KIT ↙

- Safely and rapidly lifts the sinus membrane while drilling
- Unique Stopper system that prevents over drilling into the sinus cavity
- Hydraulic Lift System that easily & safely lifts the membrane
- Bone Carrier System for transferring and filling bone material
- Bone Spreading System for spreading & compacting bone material
- Simple and intuitive surgical system
- The ability to combine Osteotome in surgery




● CAS-Drill SPECIFICATIONS & PERFORMANCE ↙


- The CAS-Drill is designed to safely and rapidly lift the maxillary sinus membrane from a crestal approach. The CAS-Drill can be used for either general-straight or tapered fixtures. It is optimized for insertion torque, initial fixation strength, and tactile feedback when using Hiossen's HG III & OSSTEM's GS / TS III Fixtures.

The CAS-Drill:


- The atraumatic design of the drill tip allows the user to perform sinus surgery even if the sinus floor is flat, incline, or septum.



- It's design forms conical bone and bone chips.
- The CAS-Drill tip has an inverse conical shape. This shape will form a conical bone chip when drilling, which assists with safely lifting the membrane. In addition, bone particles generated when drilling discharge upwards, producing a Membrane Auto-Lift function.




- Membrane can safely be lifted.




The CAS-Drill can:

- Drilling can be done at various speeds, from low to high speed (800rpm), allowing flexibility during surgery.




Guide : 400 ~ 800 rpm
However, 400 to 600rpm is recommended for first time users.

- The drill is designed with four blades which reduce deflecting off of the bone, and the straight sides dampen vibrations.



- Extraction of bone particles (at low speed of ~50rpm).




• Generally, the CAS-Drill can be used up to 50 times.

The number of uses may vary depending on the type of bone.

● Components


1) \varnothing 2.0 Twist Drill

- The drill tip is 0.6mm and is 13mm long.
- Recommended drill speed: 1000~1500 RPM
(Water Infusion + Pumping)
- 1mm spaced markers with wide bands at 4~5, 9~10
- Unique Stopper system
- It is recommended to stop drilling when there is about 2mm of bone left, please calculate this beforehand when using CT images as a guide.









2) CAS-Drill

- Comes in six (6) diameters: \varnothing 2.8 / \varnothing 3.1 / \varnothing 3.3 / \varnothing 3.6 / \varnothing 3.8 / \varnothing 4.1
- Allows a 13mm Fixture to be implanted
- Drilling is dependent upon the fixture diameter and the how far the fixture protrudes into the maxillary cavity.
- Drilling speed ranges from low speed to high speed (800rpm)
Experienced: 800rpm; Beginner: 400 to 600rpm is recommended
(Water Infusion + Pumping)
- Unique Stopper system



■ An example of a CAS-Drill dependence on the Hiossen's HG III & OSSTEM's GS / TS III Fixture diameter and protrusion height - Fixture protrusion height (mm)

Fixture	HG III, GS / TS III F \varnothing 4.0		HG III, GS / TS III F \varnothing 4.5		HG III, GS / TS III F \varnothing 5.0	
	0~3	3~6	0~3	3~6	0~3	3~6
Fixture Protrusion Height(mm)						
	\varnothing 2.8	\varnothing 3.1	\varnothing 3.3	\varnothing 3.6	\varnothing 3.8	\varnothing 4.1
CAS-Drill						
Code	SNDR2813T	SNDR3113T	SNDR3313T	SNDR3613T	SNDR3813T	SNDR4113T

● Components ↙



3) Stopper System

- A total of eleven (11) stoppers; labeled 2 to 12mm
- Labels indicate the remaining length of the drill (from drill tip to stopper top)
- Each stopper is anodized and color coded. Labels are laser etched.




4) Depth Gauge

- Measures the thickness of the remaining bone
- The atraumatic tip can be used to confirm membrane lifting
- Can be used with the Stopper system
- **Caution: Do not use the Depth Gauge to lift membrane beyond 1mm.**



5) Hydraulic Lifter

- The Hydraulic Lifter uses normal saline to raise the membrane
- **Infuse 1cc with a syringe**
- Required volume of saline
To expand 3mm of the membrane, generally 0.2 to 0.3cc of saline is injected. Inject saline very SLOWLY.
- **Contraindication**
 - Not recommended for patients with inflammation of the maxillary Sinus (Sinusitis)
 - Not recommended for patients with complex morphology of the sinus floor (including the septum)




6) Bone Carrier

- Transplanting bone material to the grafting site
- Has dual diameters: \varnothing 3.5 and \varnothing 3.9


※ Bone graft material and filler (for reference)

Herry Y and Lee DY, 2005	
Lift heigh	Volume of bone matrix
3mm	0.36cc
4mm	0.5cc
5mm	0.7cc
6mm	0.9cc



7) Bone Condenser



- Assists compacting bone grafting material
- Has dual diameters: \varnothing 2.3 and \varnothing 3.3
- Can be used with the Stopper system
- Wide banded markers at 4-5 and 9-10mm
- Can also be used to confirm membrane lifting after using the CAS-Drill



● Components ↙


8) Bone Spreader

- Evenly spreads bone material after transplanting bone material to the site.
 - After injecting 0.2 to 0.3cc use the spreader and add additional material
 - Use at lower speeds: ~ 30rpm is recommended
 - Comes in two diameters: $\varnothing 2.0$ / $\varnothing 3.0$
 - Can be used with the Stopper feature
 - The total length (head tip to stopper hilt) is 2.5mm longer other CAS-KIT tools
- Caution: When equipping this tool remember that the length is 2.5mm longer.



9) Hydraulic Lifter Tube

- Used with a syringe
- Reusable, sterilize in an autoclave.



● Clinical Indications & Case Study


1) #26, #27 Septum Case (F/36)

*Data source from: Professor Kim Gyeong-won from Chungbuk National Univ. Hospital


- #26 Septum Case
- The membrane is lifted 4~5mm and the remaining bone is about 5mm
- #26 GSII ϕ 4.0 x 10mm
- #27 GSII ϕ 4.5 x 10mm were implanted




- ϕ 2.0 Twist Drill
- Using a 3.0mm Stopper with the 2.0 Twist Drill, we are able to drill 3mm into the bone, confidently leaving 2mm of bone.



- ϕ 2.8 \rightarrow ϕ 3.1 CAS-Drill (800rpm)
- A 5mm Stopper is used for the final drilling and lifting of the membrane.



- Membrane safely lifted
- A conical bone chip is formed and pushes up the membrane, with the assistance of bone particles formed during drilling.



● Clinical Indications & Case Study

1) #26, #27 Septum Case

- Depth Gauge

Confirm membrane lifting and measuring the bone thickness



- Membrane Lift

The membrane is lifted by slowly injecting 0.30cc of saline solution using a 1cc syringe



- Bone Carrier

Osteoss Bone Powder 0.25cc is transplanted
A mix of Cortical 50%: Cancellous 50%



- Bone Condenser

Vertical compacting of the bone grafting material



- Bone Spreader

Evenly spread the bone grafting material at 10rpm of rotational speed




1) #26, #27 Septum Case


- Fixture implantation
#26 GSII \varnothing 4.0 x 10mm implanted
using 20 to 30Ncm



- Fixture implantation
#27 GSII \varnothing 4.5 x 10mm implanted
using 20 to 30Ncm



- Results
#26, #27 Missing, a case with 6mm
of bone remaining



Even though there was a Septum at #26, fixture implantation was successfully completed using the CAS-KIT to safely lift the membrane and establishing a secure implant site.

● Clinical Indications & Case Study

2) #26 Missing Case

*Data source from: Dr. So, Gwang-seup; Mirae Dental Clinic


- USII $\varnothing 4.0 \times 11.5$ mm implant planning
 - Initiated using a $\varnothing 2.0$ Twist Drill
 - CAS-Drill at 800rpm
 - Membrane lifted with 0.25cc of saline solution
 - Bone Condenser 4~5mm lifting
 - Bone Spreader at 10rpm
 - Initial fixation force 36Ncm



3) #25 Hydraulic Lift Case

*Data source from: Dr. Jung, Gi-don; Bright Smile Dental Clinic

- TSIII $\varnothing 4.5 \times 10$ mm implant planning
 - Initiated using a $\varnothing 2.0$ Twist Drill
 - CAS-Drill at 800rpm
 - Membrane lifted with 0.30cc of saline solution
 - Bone Condenser: 4mm lifting
 - Bone Spreader at 30rpm



Surgical Procedure

The CAS-Drill design is optimized for Hiossen's HG III & OSSTEM's GS / TS III Fixtures. Use the matrix below to prepare for surgery. There are a few things that need to be taken into consideration; the diameter of the fixture, the height of the fixture apex protruding into the sinus floor, and the necessary force for a stable fixture. In the case of a general straight type fixture, use a CAS-Drill that is 1mm smaller in diameter than that of the fixture.


● : Required ● : Optional

Fixture selection		Twist Drill	CAS-Drill						Hydraulic Lift & Bone Condensing					CAS-Drill
			F ϕ 4.0		F ϕ 4.5		F ϕ 5.0		Depth Gauge	Hydraulic Lifter	Bone Carrier	Bone Condenser	Bone Spreader	
Diameter(ϕ)	Fixture protruding height	ϕ 2.0	ϕ 2.8	ϕ 3.1	ϕ 3.3	ϕ 3.6	ϕ 3.8	ϕ 4.1						
F ϕ 4.0	0-3 mm	●	●						●	●	●	●	●	ϕ 3.3
	3-6 mm	●	●	●					●	●	●	●	●	-
F ϕ 4.5	0-3 mm	●	●		●				●	●	●	●	●	ϕ 3.8
	3-6 mm	●	●	●		●			●	●	●	●	●	-
F ϕ 5.0	0-3 mm	●	●		●		●		●	●	●	●	●	ϕ 4.1
	3-6 mm	●	●	●		●		●	●	●	●	●	●	-

CAS-KIT



● CAS-KIT Care & Maintenance ↙



① Prepare tools for surgery by soaking them in a “saline solution” or in “distilled water.”

② After surgery: All tools should be soaked in an “alcohol solution”.



- Avoid using Hydrogen Peroxide.
- Hydrogen Peroxide will discolor laser markings and anodized surfaces.

③ Tools should be cleaned thoroughly with distilled or tap water to wash away any remaining blood and foreign material.

④ Completely dry all tools using a dry cloth or warm air.

⑤ Dried tools should be stored in the KIT case.
(Please refer to the color coding when placing the tools back in the case)

⑥ After placing all the tools back into the kit, dry the entire kit in an Autoclave (132°C for 15 minutes) and then store the kit at room temperature.

NOTES:

It is recommended to re-sterilize the surgical KIT right before surgery. (132°C; for 15 minutes)

Immediately after surgery, all the tools should be cleaned and stored.

The CAS-KIT has a one year warranty on all parts & case.

The recommended usage of the drills is 50 times.

Crestal Approach Sinus-KIT

CAS-KIT

Edition 02/2013 EN
www.sinuskit.com
www.osstem.de

HIOSSSEN

Manufacturer
HIOSSSEN Inc.

85 Ben Fairless Dr
Fairless Hills, PA 19030, USA
TEL: 1-888-678-0001
FAX: 1-267-759-7004
www.hiossen.com

OSSTEM[®]
IMPLANT

Qualität schafft Vertrauen

Authorized Sales by
OSSTEM Germany GmbH

Mergenthalerallee 25
65760 Eschborn, Germany
Tel.: +49 (0) 6196 777 550
Fax: +49 (0) 6196 777 5529

www.osstem.de • germany@osstem.de